

## Issue 2 - Device Detection Issues

Joel Barcillo - 2024-02-09 - Troubleshooting

### Issue

No Device Found Error / Cannot sync device

### Possible Causes

- Bluetooth is not turned on
- The device was paired directly on the phone through the Bluetooth settings and not through the Welch Allyn Home app
- The device is not in "Discoverable Mode"
- The device is too far from your phone
- The initial device pairing is no longer working

### Suggested Actions

- Go to settings and make sure that Bluetooth is on and that the device is in "Discoverable" mode
- Pair the device using the Welch Allyn Home app and not under Bluetooth settings
- Please keep the device close to your mobile phone, and ensure it is turned on and in pairing mode.
- Unpair / Forget the device from the app and try the pairing process again

[Go back to the Patient Help Directory](#)

[Go back to Troubleshooting Directory](#)