



Change Units of Measure - Patient

Jay-Ileen (Ai) - 2023-07-25 - Patient Guides

[Steps on changing the units of measure](#)

You can update the units of measure in the app after initial setup by doing the following:.

1. Touch Profile.
2. Touch **Units of Measurement**.

Take note, the text below Units of measure indicates the current units selection.

1. Touch Weight, Height/Waist, and Temperature to change the units of measure.
2. Touch the desired unit of measure.
3. Touch **Save** to save the unit of measure.

[Go back to the Patient Help Directory](#)