


Customizing Supported Devices


Joel Barcillo - 2023-06-13 - Clinic Admin Guides

As a Clinic Admin, you can view and give permission to support particular Welch Allyn Home devices by following the steps below:


1. Go to the Admin Panel from your Clinic Dashboard by clicking your name and then **Admin Panel**.



2. From the Admin Panel, click on the **Clinic** tab. The Clinic tab will open showing the Clinic Details and the Clinic subtab by default.



3. From the Clinic Details screen, click on the **Health Settings** tab.



- Under the Supported Devices tab, you can see the toggle for Welch Allyn Home devices switched to **On**. Below it, You will see checkboxes for the Blood Pressure and Scale Devices with the respective names Welch Allyn Home Blood Pressure Monitor and Welch Allyn Scale 100.

CLINIC DETAILS

Clinic

Health Settings

Health Settings

Settings apply to all patients in this clinic.

4

Supported Devices

Welch Allyn Home **On**

- Blood Pressure
 - Welch Allyn Home Blood Pressure Monitor
- Scale
 - Welch Allyn Scale 100

- Warning: Please note that these settings are On and checked by default. Turning the setting off or unticking the checkboxes is not advisable as it will stop patients from being able to upload their measurements.

[Go back to the Admin Help Directory](#)