

Knowledgebase > Welch Allyn Home > Clinic Admin Guides > Health Vital Logging and Display

Health Vital Logging and Display Joel Barcillo - 2023-06-13 - Clinic Admin Guides

Further below under the Supported Devices Tab, separate boxes for each health vital are shown. Each of these items will be checked by default.

| CLINIC DETAILS | | | | | |
|---|---|----------------------------|----------------------------|----------------------------|--|
| Clinic Health Setting | 15 | | | | |
| Health Settings Settings apply to all patients | in this clinic. | | | | |
| Supported Devices | | | | | |
| Welch Allyn Home Blood Pressure Welch Allyn Home | On ^ | | | | |
| Pressure Monitor | | | | | |
| Welch Allyn Scale | 100 - | | | | |
| Velch Allyn Scale | | | | 0- | • ••••••••••••••••••••••••••••••••••• |
| Welch Allyn Scale | Blood Pressure Min. Data Entries per week 7 | Heart Rate | Min. Data Entries per week | Min. Data Entries per week | Min. Data Entries per we |
| Welch Allyn Scale | Blood Pressure Min. Data Enerles per week | Min. Data Entries per week | Min. Data Entries per week | Min. Data Entries per week | Respiratory I Min. Data Entries per we 7 Data Source |

• Warning: Please note that these settings are checked by default. Changing the settings or unticking the checkboxes is not advisable as it will stop patients from being able to upload their measurements.

Go back to the Admin Help Directory